



VitaTaq® PCR Kit & 2X Mastermix

for research use only

Procomcure Biotech GmbH
www.procomcure.com
office@procomcure.com

VitaTaq® DNA Polymerase is a robust and reliable Taq DNA polymerase suitable for all common PCR applications like colony PCR, cloning applications, high-throughput PCR and routine PCR.

PRODUCT	SIZE	SKU
VitaTaq® PCR Kit	500 units	PCCSKU1001
	1000 units	PCCSKU1002
	5000 units	PCCSKU1003
2X HS Mastermix Clear	2 ml	PCCSKU1004
	5 ml	PCCSKU1005
	10 ml	PCCSKU1006
2X HS Mastermix Gold	2 ml	PCCSKU1007
	5 ml	PCCSKU1008
	10 ml	PCCSKU1009

STORAGE CONDITIONS

Store all components at -20°C and avoid repeated freeze and thaw cycles.

ADDITIONAL MATERIALS REQUIRED

Nuclease free dH₂O

Nuclease free PCR tubes / plates & sealing options

Thermocycler

PCR Primer (10 µM each)

dNTP Mix (10 mM each)¹

template DNA

¹ (not required for 2X Mastermix Clear/Gold)

REACTION SETUP

1) Thaw all reaction components completely and mix gently to ensure even distribution of all components. Prepare the reaction on ice in a sterile, nuclease free tube and mix gently after addition of the polymerase. Collect all liquid at the bottom of the tube by a quick spin.

COMPONENT	VOLUME	FINAL CONCENTRATION
10X PCR Buffer	5 µl	1X
dNTP Mix (10 mM each)	0.5 µl	0.1 mM each
Primer 1 (10 µM)	1 µl	0.1 µM – 0.5 µM
Primer 2 (10 µM)	1 µl	0.1 µM – 0.5 µM
VitaTaq® DNA Polymerase (5 U/µL)	0.2 µl	1 U
template DNA	1 µl	<1 µg
dH ₂ O	to 50 µl	

Quality Management
System Certified

ISO 9001:2015
EN ISO 13485:2016



2) Keep the reactions on ice until transfer to the thermocycler, then cycle according to these guidelines:

STEP	CYCLES	TEMPERATURE	DURATION
Initial Denaturation	1	94°C	5 minutes
Amplification	25-35	94°C	30 seconds
		$T_m - 5^\circ\text{C}$	30 seconds
		72°C	1 minute / kb
Final Extension	1	72°C	5 minutes
Hold	1	4°C	

3) Analyze the amplification reaction by gel electrophoresis using an acrylamide or agarose gel of appropriate percentage.

2X HS Mastermix Clear / Gold

For routine applications, VitaTaq® DNA Polymerase is available as a 2x mastermix formulation; all there is to add is your primers and template DNA.

2X HS Mastermix Gold furthermore contains a tracking dye that enables direct gel loading without addition of a separate loading dye. For increased fidelity, 2X Mastermix Clear/Gold comes with a HotStart feature that allows reaction setup at room temperature.

REACTION SETUP

1) Thaw all reaction components completely and mix gently to ensure even distribution of all components. Prepare the reaction in a sterile, nuclease free tube and mix gently. Collect all liquid at the bottom of the tube by a quick spin.

COMPONENT	VOLUME	FINAL CONCENTRATION
2X HS Mastermix Clear/Gold	25 µl	1X
Primer 1 (10 µM)	1 µl	0.1 µM – 0.5 µM
Primer 2 (10 µM)	1 µl	0.1 µM – 0.5 µM
template DNA	1 µl	<1 µg
dH ₂ O	to 50 µl	

2) Transfer the reactions to the thermocycler, then cycle according to these guidelines:

STEP	CYCLES	TEMPERATURE	DURATION
Initial Denaturation	1	94°C	5 minutes
Amplification	25-35	94°C	30 seconds
		$T_m - 5^\circ\text{C}$	30 seconds
		72°C	1 minute / kb
Final Extension	1	72°C	5 minutes
Hold	1	4°C	

3) Analyze the amplification reaction by gel electrophoresis using an acrylamide or agarose gel of appropriate percentage.