



VitaProof® is a high-fidelity Pfu DNA polymerase that exhibits a 3'→5' exonuclease activity. This feature allows the removal of falsely incorporated nucleotides making VitaProof® the ideal choice for applications that require an elevated level of fidelity.

| PRODUCT | SIZE | SKU |
|--------------------|--------|------------|
| VitaProof® PCR Kit | 250 U | PCCSKU1016 |
| | 500 U | PCCSKU1017 |
| | 2500 U | PCCSKU1018 |

ADDITIONAL MATERIALS REQUIRED

- Nuclease free PCR tubes or plates
- PCR cycler
- PCR Primer
- dNTPs
- Template DNA
- Filter pipette tips
- Sterile, nuclease-free, DNA-free tubes for preparing the reaction mix

STORAGE

Store all components at -20°C and avoid repeated freeze and thaw cycles.

REACTION SETUP

- 1) Thaw all components on ice and mix gently to ensure even distribution of all components. Prepare the reaction on ice in a sterile, nuclease free tube and mix gently after addition of the polymerase. Collect all liquid at the bottom of the tube by a quick spin. Keep the reaction on ice until you transfer it to the thermocycler.

| COMPONENT | VOLUME | FINAL CONCENTRATION |
|-----------------------|--------|---------------------|
| 10X VitaProof® Buffer | 5 µl | 1X |
| dNTP Mix (10 mM each) | 1 µl | 0.2 mM each |
| Primer 1 (10 µM) | 1 µl | 0.1 µM – 0.5 µM |
| Primer 2 (10 µM) | 1 µl | 0.1 µM – 0.5 µM |
| VitaProof® 2.5 U/µl | 0.5 µl | 1.25 U |
| template DNA | 1 µl | < 1 µg |
| dH ₂ O | | to 50 µl |

- 2) Transfer the reactions to the thermocycler, then cycle according to these guidelines:

| STEP | CYCLES | TEMPERATURE | DURATION |
|-----------------|--------|----------------------|------------------|
| Initial | 1 | 94°C | 5 minutes |
| Amplification | 25-35 | 94°C | 30 seconds |
| | | T _m – 5°C | 30 seconds |
| | | 72°C | 1-2 minutes / kb |
| Final Extension | 1 | 72°C | 5 minutes |

- 3) Analyze the amplification reaction by gel electrophoresis using an acrylamide or agarose gel of appropriate percentage.