



consumables

Procomcure Biotech GmbH
Breitwies 1
5303 Thalgau - Austria
www.procomcure.com
sales@procomcure-consumables.com



10 µl



200 µl



1000 µl

Filtered Tips

10 µl, 200 µl and 1000 µl

Technical Data Sheet

Features

- filtered to prevent cross contamination
- available for 10 µl, 200 µl and 1000 µl
- non-autoclavable, tips are delivered sterile
- DNA, DNase, RNase and endotoxin-free
- racked in 10x96 tips for convenience



Pipette Compatibility

Brand	Eppendorf																		
Model	Electronic					Reference										Research		Research Pro	Research Plus
Volume	0.5-10 µl	5-100 µl	2-200 µl	50-1000 µl	0.1-2.5 µl	0.1-10 µl	0.5-10 µl	10-100 µl	20-200 µl	50-200 µl	100-1000 µl	0.1-2.5 µl	0.5-10 µl	10-100 µl	20-200 µl	100-1000 µl	20-300 µl	50-1000 µl	100-1000 µl
10 µl	•				•	•	•		•	•		•	•		•		•		
200 µl		•	•					•	•	•				•	•		•		
1000 µl				•							•					•		•	•

Brand	ThermoScientific									
Model	Finnpipette									Finnpipette Digital
Volume	0.2-2 µl	0.3-3 µl	0.2-10 µl	0.5-5 µl	0.5-10 µl	1-10 µl	5-40 µl	5-50 µl	10-100 µl	40-200 µl
10 µl	•	•	•	•	•	•	•	•	•	•
200 µl							•	•	•	•
1000 µl										•

Brand	ThermoScientific									
Model	Finnpipette electronic Multi Channel			Pipetman Classic			Pipetman Neo			Pipetman Ultra
Volume	1-10 µl	30-300 µl	100-1000 µl	P2, 0.1-2 µl	P2, 0.2-2 µl	P10, 0.2-10 µl	P20, 2-20 µl	P100, 10-100 µl	P200, 20-200 µl	P1000, 100-1000 µl
10 µl	•	•	•	•	•	•	•	•	•	•
200 µl										
1000 µl			•							

Brand	VWR							
Model	VWR & VWR Signature							
Volume	0.1-2 µl	0.5-10 µl	2-20 µl	5-50 µl	10-100 µl	20-200 µl	50-250 µl	100-1000 µl
10 µl	●	●						
200 µl			●	●	●	●	●	
1000 µl								●

Brand	Starlab		
Model	ErgoOne		
Volume	0.5-10 µl	10-100 µl	100-1000 µl
10 µl	●		
200 µl		●	
1000 µl			●

Brand	Axygen						
Model	AxyPet						
Volume	0.1-2 µl	0.5-10 µl	2-20 µl	5-50 µl	10-100 µl	20-200 µl	100-1000 µl
10 µl	●	●					
200 µl			●	●		●	
1000 µl							●

Brand	BrandTech						
Model	Transferpette						
Volume							
	2-20 µl	0.5-10 µl	5-50 µl	10-100 µl	20-200 µl	25-250 µl	100-1000 µl
10 µl		●					
200 µl	●		●	●	●	●	
1000 µl							●

Brand	Biohit									
Model	Proline					m-Line			e-Line	
Volume	0.1-2.5 µl	0.2-10 µl	0.5-10 µl	5-50 µl	5-100 µl	10-250 µl	10-500 µl	20-200 µl	50-1000 µl	100-1000 µl
10 µl	•	•	•	•	•	•	•	•	•	•
200 µl										
1000 µl										

Brand	Rainin									
Model	EDP3					Latch Mode				
Volume	E-MIC-10, 0.5-10 µl	E2-25, 2.5-25 µl	E2-100, 10-100 µl	E2-250, 25-250 µl	E2-1000, 100-1000 µl	R2, 0.1-2 µl	R10, 0.2-10 µl	R20, 2-20 µl	R100, 10-100 µl	R200, 20-200 µl
10 µl	•	•	•	•	•	•	•	•	•	•
200 µl										
1000 µl										

Contact & Ordering

Filtered Tips, 10 µl, racked & sterile
 Filtered Tips, 200 µl, racked & sterile
 Filtered Tips, 1000 µl, racked & sterile

10 x 96 tips
 10 x 96 tips
 10 x 96 tips

PCCSKU30003
 PCCSKU30004
 PCCSKU30005

For ordering and inquiries, contact us at sales@procomcure-consumables.com



www.procomcure.com