



dNTP Mix (10 mM each)

for research use only

Procomcure Biotech GmbH
www.procomcure.com
office@procomcure.com

dNTP Mix (10 mM each)

PRODUCT	SIZE	SKU
dNTP Mix (10 mM each)	0.5 ml	PCCSKU1019
	1 ml	PCCSKU1020
	2.5 ml	PCCSKU1021

STORAGE & HANDLING

Store at -20°C. Thaw and mix properly before use. Avoid excessive freeze and thawing cycles, prepare aliquots if necessary.

DESCRIPTION

dNTPs are an essential component of every PCR reaction. dATP, dTTP, dCTP and dGTP are supplied in a convenient ready-to-solution with 10 mM of each nucleotide that can be used directly for most PCR protocols without further dilution.

FEATURES

- Ultrapure: $\geq 99\%$ as determined by HPLC
- Supplied in ultrapure dH₂O
- Nuclease-free
- Highly stable
- Compatible with all Procomcure PCR reagents e.g. VitaTaq® DNA Polymerase